

Akkommodative® 1CU



Technical Information


● Accommodative Foldable Monobloc Posterior Chamber Lens

● Material:

- Hydrophilic MicroCryl®
- UV inhibitor


● Est. A-Constant: 118.1*

(*for current status please refer to www.humanoptics.com)



Diopter Range

biconvex



In 0.5 D steps

16.0 D to 35.0 D

Parameters for Lens Calculation with AkkoWizard

Target refraction

-0.5 D to 0.0 D

Calculation formula

Haigis

Optimised constants

a0: 1.833 / a1: 0.225 / a2: 0.094

Indication Recommendations

- ✓ Cataract
- + Enhanced Reading Comfort



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by