

VR 23 G

by Janach



Janach

J 3834.23L

*ILM Internal Limitant Membrane forceps,
23 G, new handle with colour code, Ø 7.5 mm
(based on Eckardt design)*



Janach

J 3839.23L

*Retinal forceps,
23 G, straight, serrated jaws,
new handle with colour code, Ø 7.5 mm*



Janach

J 3870.23L

*Retinal scissors,
23 G, slightly curved jaws, new handle with colour code, Ø 7.5 mm
To be used with Trocar*

*Strumenti chirurgici per oftalmologia
e microchirurgia*

Ophthalmic Surgical Instruments

CE

*Sistema di Qualità
Certificato*

*Certified
Quality System*

*ISO 9001
ISO 13485*

*© Copyright:
Hermann Janach*

*E. Janach srl
Via Borgo Vico, 35 - 22100 Como - Italy
Tel. +39 031 574088 - Fax +39 031 572055
international@janach.it - www.janach.it*

e.janach®



VR EVOLUTION

by Janach


Janach
J 3009.10

*Sideport incisions multipurpose microscissors,
23 G, for anterior chamber,
retinal surgery,
titanium handle*




Janach
J 3867

*18 G foreign body forceps direct action,
with 5 mm tips,
tungsten carbide coated tips,
for secure gripping*



Janach
J 3868

*Magnet microforceps for foreign body,
positive action*



*Strumenti chirurgici per oftalmologia
e microchirurgia*

Ophthalmic Surgical Instruments

CE

*Sistema di Qualità
Certificato*

*Certified
Quality System*

*ISO 9001
ISO 13485*

© Copyright:
Hermann Janach

*E. Janach srl
Via Borgo Vico, 35 - 22100 Como - Italy
Tel. +39 031 574088 - Fax +39 031 572055
international@janach.it - www.janach.it*

e.janach®

