

Aspira[®] -aA



Technical Information

- **Foldable Monobloc Posterior Chamber Lens**

- **Material:**


- Hydrophilic MicroCryl[®]
- UV inhibitor

- **Optic:**

- Aberration free (aspheric posterior surface)
- 360° lens epithelial cell barrier


- **Est. A-Constant:** 118.4*

(*for current status please refer to www.humanoptics.com)



Diopter Range

biconvex




In 0.5 D steps

10.0 D to 30.0 D

In 1.0 D steps

0.0 D to 9.0 D

concave/convex



Indication Recommendations

✓ Cataract



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by