

Aspira[®] -aAY



Technical Information

● Foldable Monobloc Posterior Chamber Lens

● Material:


- Hydrophilic MicroCryl[®]
- UV inhibitor
- Yellow colouring agent for blue light protection

● Optic:

- Aberration free (aspheric posterior surface)
- 360° lens epithelial cell barrier


● Est. A-Constant: 118.4*

(*for current status please refer to www.humanoptics.com)



Diopter Range

biconvex




In 0.5 D steps

10.0 D to 30.0 D

In 1.0 D steps

0.0 D to 9.0 D

concave/convex



Indication Recommendations

- ✓ Cataract
- + Retina Protection



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by