

DiffRACTIVA® Diff-sPBY




Technical Information

• Foldable Diffractive 3-Piece Posterior Chamber Add-On Lens for Sulcus Fixation in Pseudophakic Eyes

• Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection
- Concave posterior surface
- Round anterior optic edge to prevent iris irritation




• Near Addition: + 3.5 D (IOL plane)

• Haptic:

High Molecular PMMA, enlarged C-Loop

Diopter Range

convex / concave



Sphere

0.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations



- ✓ Pseudophakic Ametropia
- + Enhanced Reading Comfort
- + Retina Protection

Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen, Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by