

DiffRACTIVA® Diff-sPB



Technical Information

• Foldable Diffractive 3-Piece Posterior Chamber Add-On Lens for Sulcus Fixation in Pseudophakic Eyes


• Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Concave posterior surface
- Round anterior optic edge to prevent iris irritation

• Near Addition: + 3.5 D (IOL plane)


• Haptic:

High Molecular PMMA, enlarged C-Loop



Diopter Range

convex / concave



Sphere

0.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Pseudophakic Ametropia
+ Enhanced Reading Comfort



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen, Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by