

Diffactiva® Diff-s



Technical Information

- **Foldable Diffractive 3-Piece Posterior Chamber Lens**

- **Optic:**

- Hydrophobic MicroSil®
- UV inhibitor
- 360° sharp optic edge


- **Near Addition: +3.5 D**

- **Haptic:**

MicroPlex®, modified C-Loop


- **Est. A-Constant: 118.6***

(*for current status please refer to www.humanoptics.com)



Diopter Range

biconvex



In 0.5 D steps

10.0 D to 30.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ **Cataract**
+ Enhanced Reading Comfort



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by