

Diffactiva® Diff-SS



Technical Information

● Foldable Diffractive 3-Piece Posterior Chamber Lens for Sulcus Fixation

● Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Round anterior optic edge to prevent iris irritation


● Near Addition: + 3.5 D

● Haptic:

MicroPlex®, enlarged C-Loop


● Est. A-Constant: 118.0*

(*for current status please refer to www.humanoptics.com)



Dioppter Range

biconvex



In 0.5 D steps

available
on request

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Aphakia or Trauma
- + Enhanced Reading Comfort



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by