



## Technical Information

- **Foldable 3-Piece Posterior Chamber Lens**

- **Optic:**


- Hydrophobic MicroSil®
- UV inhibitor
- Aberration correcting (aspheric posterior surface)
- 360° sharp optic edge

- **Haptic:**

MicroPlex®, modified C-Loop


- **Est. A-Constant:** 118.6\*

(\*for current status please refer to [www.humanoptics.com](http://www.humanoptics.com))



## Diopter Range

### biconvex



### In 0.5 D steps


7.0 D to 31.0 D

### In 1.0 D steps

32.0 D to 36.0 D

-6.0 D to 6.0 D

### concave/convex



## Customized Lenses

Other diopters can be delivered according to individual specifications

## Indication Recommendations



Cataract

+ Extreme Hyperopia

+ Extreme Myopia

+ Improved Functional Vision

+ Reduced Glare Sensitivity

### Manufactured by

HumanOptics AG  
Spardorfer Str. 150  
91054 Erlangen · Germany  
Tel: + 49 91 31/50 66 5.0  
Fax: + 49 91 31/50 66 5.90  
[mail@humanoptics.com](mailto:mail@humanoptics.com)  
[www.humanoptics.com](http://www.humanoptics.com)

### Distributed by