

Secūra® -sSY



Technical Information

- **Foldable 3-Piece Posterior Chamber Lens for Sulcus Fixation**

- **Optic:**


- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection

- **Haptic:**

MicroPlex®, enlarged C-Loop


- **Est. A-Constant:** 118.0*

(*for current status please refer to www.humanoptics.com)



Dioppter Range

biconvex



In 0.5 D steps

7.0 D to 31.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Aphakia or Trauma
+ Retina Protection



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by