




## Technical Information

### • Foldable 3-Piece Posterior Chamber Add-On Lens for Sulcus Fixation in Pseudophakic Eyes

#### • Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Concave posterior surface
- Yellow colouring agent for blue light protection
- Round anterior optic edge to prevent iris irritation




#### • Haptic:

- High Molecular PMMA
- Enlarged C-Loop

## Diopter Range

convex / concave



In 0.5 D steps

-6.0 D to 6.0 D

## Customized Lenses

Other diopters can be delivered according to individual specifications

## Indication Recommendations

- ✓ Pseudophakic Ametropia + Retina Protection



### Manufactured by

HumanOptics AG  
Spardorfer Str. 150  
91054 Erlangen, Germany  
Tel: + 49 91 31/50 66 5.0  
Fax: + 49 91 31/50 66 5.90  
mail@humanoptics.com  
www.humanoptics.com

### Distributed by