

Torica® -sPBV




Technical Information

• Foldable Toric 3-Piece Posterior Chamber Add-On Lens for Sulcus Fixation in Pseudophakic Eyes

• Optic:


- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection
- Toric anterior surface
- Concave posterior surface



• Haptic:

High Molecular PMMA, undulated, enlarged C-Loop

Diopter Range



Sphere

Adjusted to cylindrical power so that spherical equivalent = 0D

Cylinder

In 1.0 D steps

1.0 D to 6.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Pseudophakic Ametropia
- + Corneal Astigmatism
- + Retina Protection



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen, Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by