

Torica® -sPB



Technical Information


• Foldable Toric 3-Piece Posterior Chamber PiggyBack Lens for Sulcus Fixation

• Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Toric anterior and spheric, concave posterior surface


• Haptic:

MicroPlex®, undulated, enlarged C-Loop



in mm

Diopter Range



Cylinder

In 1.0 D steps

1.0 D to 6.0 D

Sphere

Adjusted to cylindrical power so that spherical equivalent = 0D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Pseudophakic Ametropia + Corneal Astigmatism



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by