




## Technical Information

- **Foldable Toric 3-Piece Posterior Chamber PiggyBack Lens for Sulcus Fixation**

- **Optic:**


- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection
- Toric anterior and spheric, concave posterior surface



- **Haptic:**

MicroPlex®, undulated, enlarged C-Loop

## Dioptr Range



### Sphere

Adjusted to cylindrical power so that spherical equivalent = 0D

### Cylinder

In 1.0 D steps

1.0 D to 6.0 D

## Customized Lenses

Other diopters can be delivered according to individual specifications

## Indication Recommendations



- ✓ Pseudophakic Ametropia
- + Corneal Astigmatism
- + Retina Protection

### Manufactured by

HumanOptics AG  
Spardorfer Str. 150  
91054 Erlangen · Germany  
Tel: + 49 91 31/50 66 5.0  
Fax: + 49 91 31/50 66 5.90  
mail@humanoptics.com  
www.humanoptics.com

### Distributed by