



Technical Information

- **Foldable Toric 3-Piece Posterior Chamber Lens for Sulcus Fixation**

- **Optic:**


- Hydrophobic MicroSil®
- UV inhibitor
- Spheric anterior and toric posterior surface

- **Haptic:**

MicroPlex®, enlarged C-Loop


- **Est. A-Constant: 118.0***

(*for current status please refer to www.humanoptics.com)



in mm

Dioppter Range



anterior

in mm

Sphere

In 0.5 D steps

15.0 D to 25.0 D

In 1.0 D steps

-3.0 D to 14.0 D

26.0 D to 31.0 D

Cylinder

In 1.0 D steps

2.0 D to 12.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Aphakia or Trauma + Corneal Astigmatism



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by