

Torica® -SSY



Technical Information

● Foldable Toric 3-Piece Posterior Chamber Lens for Sulcus Fixation

● Optic:


- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection
- Spheric anterior and toric posterior surface

● Haptic:


MicroPlex®, enlarged C-Loop

● Est. A-Constant: 118.0*

(*for current status please refer to www.humanoptics.com)



Diopter Range



Sphere

In 0.5 D steps

15.0 D to 25.0 D

In 1.0 D steps

-3.0 D to 14.0 D

26.0 D to 31.0 D

Cylinder

In 1.0 D steps

2.0 D to 12.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations

- ✓ Aphakia or Trauma
- + Corneal Astigmatism
- + Retina Protection



Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by