

Torica® -SY




Technical Information

● Foldable Toric 3-Piece Posterior Chamber Lens

● Optic:

- Hydrophobic MicroSil®
- UV inhibitor
- Yellow colouring agent for blue light protection
- Spheric anterior and toric posterior surface
- 360° sharp optic edge




● Haptic:

MicroPlex®, undulated Z-Loop

● Est. A-Constant (Haigis): 118.6*

(*for current status please refer to www.humanoptics.com)

Diopter Range



Sphere

In 0.5 D steps

15.0 D to 25.0 D

In 1.0 D steps

-3.0 D to 14.0 D

26.0 D to 31.0 D

Cylinder

In 1.0 D steps

2.0 D to 12.0 D

Customized Lenses

Other diopters can be delivered according to individual specifications

Indication Recommendations



- ✓ Cataract
- + Corneal Astigmatism
- + Extreme Hyperopia
- + Extreme Myopia
- + Retina Protection

Manufactured by

HumanOptics AG
Spardorfer Str. 150
91054 Erlangen · Germany
Tel: + 49 91 31/50 66 5.0
Fax: + 49 91 31/50 66 5.90
mail@humanoptics.com
www.humanoptics.com

Distributed by